en Pre-lubricated urethral intermittent catheter with integrated collection bag

Contents

Pre-lubricated urethral intermittent catheter with integrated collection bag

Product Description:

- Ready-to-use, pre-lubricated with gel, single use intermittent catheter
- Protective tip and integrated urine collection bag support no touch insertion technique
- Gel reservoir
- 2 smooth catheter eyelets
- Colour-coded funnel

	Catheter effective length: the length of the catheter that can be inserted into the urethra.		
40 cm/ 16" catheter	39 cm	390 mm	15.4"
20 cm/ 8" catheter	19 cm	190 mm	7.5"

Not made with natural rubber latex. This product consists of: polyethylene/EVA collection bag, thermoplastic elastomer tip, water soluble lubricating gel, and tubing not made with DEHP.
Indications and Intended Use: This intermittent catheter is a flexible tubular device that is inserted through the urethra by male, female and paediatric patients who need to drain urine from the bladder.

WARNINGS

Do not use if package is damaged or opened. To help reduce the potential for infection and/or other complications, do not reuse. If discomfort or any sign of trauma occurs, discontinue use immediately and consult your healthcare practitioner.


Precautions: Please consult your healthcare practitioner before using this product if any of the following conditions are present:

- Severed urethra
- Unexplained urethral bleeding
- Pronounced stricture
- False passage
- Urethritis – inflammation of the urethra
- Prostatitis – inflammation of the prostate gland
- Epididymitis – inflammation of the epididymis (testicle tube)

Self-catheterisation should follow the plan of care and advice given by your healthcare practitioner and be carried out only in accordance with the instructions for use provided. Because catheterisation frequency varies by person, the recommended frequency of your catheterisation should be provided by your healthcare practitioner. For any other questions about your catheterisation, please contact your healthcare practitioner.

NOTES: Contact a healthcare practitioner for help with continence care products. In case of serious injury (incident) in relation to your use of the product, please contact your local distributor or manufacturer, and your local competent authority. For more information, see www.hollister.com/authority or contact EC Rep or local distributor.

Please see instructions section.

			
Urethral Intermittent catheter			
	Single sterile barrier system with protective packaging outside		
	Single sterile barrier system		
	Consult Instructions for Use if package is damaged		Number of units
	CE mark as required		Medical Device (as applicable)

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