	en 2 Piece Drainable Pouch
	da Tømlar 2-dels pose
	de Zweiteiliger Ausstreifbeutel
	es Bolsa abierta de 2 piezas
	fi 2-osainen tyhjennettävä avanpessu
	fr Poche visible 2 pièces
	fr-CA Sac visible 2 pièces
	it Sacca per ileostomia per sistema due pezzi
	nl Tweedelig open opvangzakje
	no Todelt tømlar pose
	sv 2-dels tømlar påse

www.hollister.com/symbolglossary

912756

Hollister Incorporated
2000 Hollister Drive
Libertyville, IL 60040 USA

Hollister LLC
Foxford Road, Ballina
Co. Mayo, Ireland

©2020 Hollister Incorporated.
912756-220

en INSTRUCTIONS FOR USE

Contents

- Drainable pouch with attached flange
- Ostomy belt tabs (belt not included)
- Integrated Lock 'n Roll™ microseal closure or clamp closure (see product label)
- Disposal bags

Product Description: Part of a two-piece ostomy pouching system. Designed to be used with a corresponding skin barrier with the same flange size. The pouch is designed to be emptied. The pouch is most suitable for colostomies or ileostomies. Product is made of polymer/copolymer plastics, and nonwoven fabric.

Indications and Intended Use: This product is part of a two-piece ostomy pouching system intended to collect output from a stoma (e.g., colostomy, ileostomy) or skin opening (e.g., wound, fistula).

WARNINGS

Reuse by reprocessing, cleaning, disinfection, and sterilisation (sterilization) may compromise the product. This may result in odour (odor) or leakage.

Notes: Contact a stoma care nurse or healthcare practitioner for help with the ostomy pouching system or skin irritation. In case of serious injury (incident) in relation to your use of the product, please contact your local distributor or manufacturer, and your local competent authority. For more information, see www.hollister.com/authority or contact EC Rep or local distributor.

MD Medical Device
(as applicable).

Not made with natural rubber latex.
CE mark as required

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